

MEDICAL CARE PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS





Empower Your Mobility Discover Elective and Manual Wheelchairs









www.dminteract.com

info@dminteract.com

801 Dundas Street East, Unit #D6, Mississauga 📀







PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



DM-MWC-01

DMInteract® DM-MWC-01

Handicap Manual Wheelchair

The DMInteract **DM-MWC-01**Foldable Lightweight Handicapped Manual Wheelchair is designed to deliver a superior mobility experience, combining strength, comfort, and ease of use.





- Frame Material: Steel Pipes
- Tyres: 6" Solid Front Wheel & 16" Rear Wheel
- Max Load: 100 Kgs
- Folded Width: 23 cm
- Seat Width: 44 cm
- · Seat Height: 47 cm
- Seat Depth: 39 cm
- Backrest Height: 42 cm
- Weight: 12.00 Kg
- Dimensions: 78 x 58 x 84 cm
- Package Weight: 13 Kgs
- Packaging Dimensions: 79 x 26 x 72 cm



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People

DM-MWC-02



Handicap Manual Wheelchair

The DMInteract **DM-MWC-02**Foldable Lightweight Handicapped Manual Wheelchair is designed to deliver a superior mobility experience, combining strength, comfort, and ease of use.





- Drive Mechanism: Manual Wheelchair
- Frame Material: Steel Pipes
- Tyres: 6" Front Wheel & 20" Rear Wheel
- Max Load: 100 Kgs
- Folded Width: 28 cm
- Seat Width: 45 cm
- Seat Height: 47 cm
- Seat Depth: 39 cm
- Backrest Height: 42 cm
- Weight: 12.00 Kg
- Dimensions: 78 x 63 x 84 cm
- Package Weight: 13 Kgs





PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



DM-MWC-03

DMInteract® DM-MWC-03

Handicap Manual Wheelchair

The DMInteract **DM-MWC-03**Foldable Lightweight Handicapped Manual Wheelchair combines functionality, comfort, and style to offer a superior mobility solution.



- Drive Mechanism: Manual Wheelchair
- Frame Material: Steel Pipes
- Tyres: 6" Front Wheel & 20" Rear Wheel
- Max Load: 100 Kgs
- Folded Width: 28 cm
- Seat Width: 45 cm
- Seat Height: 47 cm
- Seat Depth: 39 cm
- Backrest Height: 42 cm
- Weight: 12.00 Kg
- Dimensions: 78 x 63 x 84 cm
- Package Weight: 13 Kgs



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



ROAD BUSTER D1

DMInteract® OPTIMUS D1

Handicap Electric Wheelchair

The DMInteract **Optimus D1** Foldable Lightweight Handicapped Electric Wheelchair is designed to deliver a superior mobility experience, combining strength, comfort, and ease of use.





- Front wheel: 6-inch solid front wheel.
- Rear wheel: 16-inch solid rear wheel.
- Seat cushion: breathable bee net cushion.
- Footrest: Foldable footrest
- Load: 100 kg
- Dimensions: width: 58cm Hight: 23cm
- Seat Width: 44cm
- Sitting Height: 47cm
- Vehicle Height: 84cm
- Vehicle Length: 78cm
- Seat Depth: 39cm
- Backrest Height: 42cm
- G.W: 13KG



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



OPTIMUS D1

DMInteract® OPTIMUS D1

Handicap Electric Wheelchair

The DMInteract **Optimus D1** Foldable Lightweight Handicapped Electric Wheelchair is designed to deliver a superior mobility experience, combining strength, comfort, and ease of use.





- Frame Material: Steel Pipes
- Motor Power: 24V/250W x 2pcs Brush Motor
- Battery Type: Lead-Acid 24v12.8Ah
- Tyres: 10" & 16" PU or Pneumatic Tire
- Max Load: 100 Kgs
- Speed: 6KM/H
- Driving Range: 15-20KM
- Folded Width: 35.5 cm
- Seat Width & Height: 45 cm x 44 cm
- Seat Depth: 46 cm
- Backrest Height: 44 cm
- Weight: 45.00 Kg
- Dimensions: 108.5 x 68.5 x 91 cm



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



OPTIMUS D2

DMInteract® OPTIMUS D2

Handicap Electric Wheelchair

The DMInteract **Optimus D2**Foldable Lightweight Handicapped Electric Wheelchair is designed to deliver a superior mobility experience, combining strength, comfort, and ease of use.





- Frame Material: Steel Pipes
- Motor Power: 24V/250W x 2pcs Brush Motor
- Battery Type: Lead-Acid 24v12.8Ah
- Tyres: 10" & 16" PU or Pneumatic Tire
- Max Load: 100 Kgs
- Speed: 6KM/H
- Driving Range: 15-20KM
- Folded Width: 35.5 cm
- Seat Width & Height: 45 cm x 44 cm
- Seat Depth: 46 cm
- Backrest Height: 44 cm
- Weight: 45.00 Kg
- Dimensions: 108.5 x 68.5 x 91 cm



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



BUMBLEBEE D1

DMInteract® BUMBLEBEE D1

Handicap Electric Wheelchair

Introducing the DMInteract

Bumblebee D1 Foldable Lightweight
Electric Wheelchair, expertly crafted
to enhance mobility and
independence for individuals with
mobility challenges.



- Frame Material: Steel Pipes
- Motor Power: 24V/250W x 2pcs Brush Motor
- Battery Type: Lead-Acid 24v12.8Ah
- Tyres: 10" & 24" PU or Pneumatic Tire
- Max Load: 100 Kgs
- Speed: 6KM/H
- Driving Range: 15-20KM
- Folded Width: 35.5 cm
- Seat Width & Height: 45 cm 44 cm
- Seat Depth: 46 cm
- Backrest Height: 44 cm
- Product Weight: 45 Kg
- Dimensions: 117.5 x 68.5 x 91 cm
- Packaging Weight: 46.90 KG



PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



BUMBLEBEE D2

DMInteract® BUMBLEBEE D2

Handicap Electric Wheelchair

The DMInteract Bumblebee D2
Foldable Lightweight Electric
Wheelchair is a perfect blend of comfort, functionality, and modern design for those seeking independence and mobility.





KEY FEATURES

• Frame: Steel Pipes

Motor Power :24V/250W*2pcs Brush Motor

Battery: Lead-acid 24v12.8Ah

• Tyres: 10" & 16" PU or Pneumatic Tire

Max Load: 100kgs

• Speed: 6KM/H

• Range: 15-20KM

• Overall Width & Height: 68cm x 89cm

Folded Width: 35.5cm

Seat Width: 45cm

Seat Height: 44cm

• Seat Deepth: 46cm Backrest Height: 44cm



FROM VISULIZATION TO COLLABORATION

MEDICAL CARE

PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



BUMBLEBEE D3

DMInteract® BUMBLEBEE D3

Handicap Electric Wheelchair

The DMInteract **Bumblebee D3**Foldable Lightweight Handicapped
Electric Wheelchair is a perfect
combination of comfort, power, and
practicality for individuals with
mobility needs.



- Frame: Steel Pipes
- Motor Power :24V/250W*2pcs Brush Motor
- Battery: Lead-acid 24v12.8Ah
- Tyres: 10" & 24" PU or Pneumatic Tire
- Max Load: 100kgs
- Speed: 6KM/H
- Range: 15-20KM
- Overall Width & Height: 68cm x 89cm
- Folded Width: 35.5cm
- Seat Width: 45cm
- Seat Height: 44cm
- Seat Deepth: 46cm Backrest Height: 44cm



FROM VISULIZATION TO COLLABORAT

MEDICAL CARE PRODUCT CATALOGUE

MOBILITY & ACCESSIBILITY AIDS

Used for Patients with Inconvenient Mobility or Handicapped People



BUMBLEBEE D4

DMInteract® BUMBLEBEE D4

Handicap Electric Wheelchair

The DMInteract **Bumblebee D4**Foldable Lightweight Handicapped
Electric Wheelchair is engineered to
offer maximum mobility and comfort
for individuals with physical
challenges.



- Frame: Steel Pipes
- Motor Power :24V/250W*2pcs Brush Motor
- Battery: Lead-acid 24v12.8Ah
- Tyres: 10" & 16" / 24" PU or Pneumatic Tire
- Max Load: 100kgs
- Speed: 6KM/H
- Range: 15-20KM
- Overall Width & Height: 45cm x 44cm
- Folded Width: 35.5cm
- Seat Width: 45cm
- Seat Height: 44cm
- Seat Deepth: 46cm Backrest Height: 44cm